

# AFLP Protocol

27 June 2014

Anthony J Geneva

Modified from:

Vos, P., R. Hogers, M. Bleeker, M. Reijans, T. Lee, M. Hornes, et al. 1995. AFLP: a new technique for DNA fingerprinting. *Nucleic Acids Research* 23:4407-4414.

and

<https://www.msu.edu/user/hazensam/aflp/AFLPprotocolMSU.html>

## Oligos

### Adaptors

EcoR1adapterA CTCGTAGACTGCGTACC  
EcoR1adapterB AATTGGTACGCAGTCTAC  
Mse1adapterA GACGATGAGTCCTGAG  
Mse1adapterB TACTCAGGACTCAT

### Pre Selective Primers

Prese1E GACTGCGTACCAATTCA  
Prese1M GATGAGTCCTGAGTAAC

### Mse Selective Unlabeled Primer

SelampM47 GATGAGTCCTGAGTAACAA  
SelampM51 GATGAGTCCTGAGTAACCA  
SelampM55 GATGAGTCCTGAGTAACGA  
SelampM59 GATGAGTCCTGAGTAACTA  
SelampM49 GATGAGTCCTGAGTAACAG  
SelampM53 GATGAGTCCTGAGTAACCG  
SelampM61 GATGAGTCCTGAGTAACTG

### EcoR1 Selective 5'-Labeled Primer

SelAmpE1 GACTGCGTACCAATTCAGG  
SelAmpE2 GACTGCGTACCAATTCATC  
SelAmpE3 GACTGCGTACCAATTCACA

Dilute all stocks to 100 $\mu$ M and store -20°C.

## Adaptor Construction

Component	Amount (ml)
Mse1adapterA (100 $\mu$ M)	7.5
Mse1adapterB (100 $\mu$ M)	7.5
dH <sub>2</sub> O	135
<b>Total (5<math>\mu</math>M)</b>	<b>150</b>

Component	Amount (ml)
EcoR1adapterA (100 $\mu$ M)	75
EcoR1adapterB (100 $\mu$ M)	75
<b>Total (50<math>\mu</math>M)</b>	<b>150</b>

### Cycler Protocol

95°C 5 mins

**70X** 94°C 30 secs/ -1°C per cycle

4°C hold

Store at -20C

### Digestion

Component	Amount ( $\mu$ l)	Per plate	/2
dH <sub>2</sub> O	13.6	1428	714
NEB 10X CutSmart	2	210	105
NEB EcoR1 (20 units/ $\mu$ l)	0.25	26.25	13.125
NEB Mse1 (10 units/ $\mu$ l)	0.05	5.25	2.625
NEB 100X BSA (20mg/mL)	0.1	10.5	5.25
DNA (50ng/ $\mu$ l)	4		
<b>Total</b>	<b>20</b>		

### Cycler Protocol

37°C 180 mins

60°C 15 mins

### Ligation

Component	Amount ( $\mu$ l)	Per plate
dH <sub>2</sub> O	5.8	609
10X T4 ligase buffer (with ATP)	2	210
EcoR1 adapter (5 $\mu$ M)	1	105
Mse1 adapter (50 $\mu$ M)	1	105
T4 DNA ligase (400 units/ $\mu$ l)	0.2	21
Entire Digestion Product	20	
<b>Total</b>	<b>30</b>	

Incubate overnight at 37°C.

Run 4 $\mu$ L of Ligation product on a 1.5% agarose gel. A diffuse smear (or sometimes distinct bands) should be visible between ~200-1000 bp. If these appear, perform a full of plate preselective amplifications. If not, test a small number of samples via preselective amplification and continue if these are successful.

### Preselective Amplification

Component	Amount ( $\mu$ l)	Per plate	/4
dH <sub>2</sub> O	22.8	2394	598.5
E Primer (10 $\mu$ M)	5	525	131.25
M Primer (10 $\mu$ M)	5	525	131.25
MgSO <sub>4</sub>	5	525	131.25
10x BioBasic Buffer	5	525	131.25
dNTPs (0.5mM)	5	525	131.25
Taq	0.25	26.25	6.5625
Ligation Product	2		
<b>Total</b>	<b>50.05</b>		

### Cycler Protocol

94°C 120 secs

94°C 60 secs

**26x** 56°C 60 secs

72°C 60 secs

72°C 60 secs

4°C hold

Run 4 $\mu$ L of Preselective Product on a 1.5% agarose gel. A DNA smear or distinct bands should appear between 50-500 bp. Store product at -20°C.

### Selective Amplification

Component	Amount ( $\mu$ l)	Per Plate	/2
H <sub>2</sub> O	11.4	1197	598.5
10x BioBasic Buffer	2.5	262.5	131.25
dNTPs (0.5mM)	2.5	262.5	131.25
Mse1 selective primer (2 $\mu$ M)	2.5	262.5	131.25
EcoR1 labeled selective primer (2 $\mu$ M)	2.5	262.5	131.25
MGSO <sub>4</sub>	2.5	262.5	131.25
Taq	0.125	13.125	6.5625
Template	1		
<b>Total</b>	<b>25.025</b>		

### Cycler Protocol

94°C 120 secs

94°C 30 secs

**12x** 65°C 30 secs/ -1°C per cycle

72°C 60 secs

94°C 30 secs

**23x** 56°C 30 secs

72°C 60 secs

72°C 60 secs

4°C hold

Run 4 $\mu$ L of Selective Product. Bands should appear between 50-500 bp. Store product at 4°C in the dark by wrapping in foil.

## Analysis Pipeline

Submit Plate to Genomics Center

### Download Data from FGC

- Genotyping tab, click on Download Your Fluorescent DNA Fragment Analysis Results
  - Unzip files

### Process data in PeakScanner

- export data

### Further process in RawGeno

- export data
- Use conversion script 'RawGeno\_2\_AFLPScore'

### AFLPScore

- Determine locus and phenotype thresholds
  - Calculate error rates
  - Export final matrix